

NJ Phillips Plastic Semi-Automatic Injector

50ml (PAS803)(PAS1357 polyprop barrel)

25ml (PAS802)(PAS1356 polyprop barrel)



OPTIONAL POLYPROPYLENE BARREL

*The instruments are equipped as standard with a very high strength, crystal clear, long lasting barrel which is suitable for use with most animal health formulations. If you experience any chemical damage to the standard barrel, we recommend the optional barrel kit WX1719 (50ml) or WX1718 (25ml) made from polypropylene.

MODEL SHOWN: PAS1357



EASY PRODUCT IDENTIFICATION
INTERCHANGEABLE COLOURED PUSH KNOB (WX1421)

MODEL SHOWN: PAS802

TIME SAVER
BUILT IN NEEDLE MAGAZINE

Suitable for: Cattle, Pigs & Sheep

- Easy choice of 5 fixed dose settings.
- Cost effective, reliable and durable to save you money.
- Built in needle magazine saves you time.
- Ergonomic fingertip action reduces hand strain.
- PAS1357 & PAS1356 come equipped with a polypropylene barrel.

Spare Parts Kits:

Available through our distributors.

- WX1421 – Pack of 4 interchangeable coloured push knobs for injectable product identification.
- WX1433 – Pack of 5 O-rings (50ml).
- WX1434 – Pack of 5 O-rings (25ml).
- Cylinder Assembly Kits (see below).

LEGEND:

1. Luer Lock Needle Nut (SN402)	15. Right Handle (SH880R/G4)
2. Needle Nut Seal Ring (SP69)	16. Needle Magazine (SH885)
3. Fibre Seal Washer (SP42)	17. Ratchet Plate Spring (SP452)
4. 50ml Cylinder (SH393)	18. Ratchet Plate (SN247)
5. 25ml Cylinder (SH823)	19. 50ml Dose Adjustor (SH437)
6. Delivery Cage Seal Ring (SP69)	20. 25ml Dose Adjustor (SH822)
7. Luer Lock Mount Delivery Cage (SN401)	21. Dose Adjustor Locking Sleeve (SH436)
8. Push Rod Guide (SH435)	22. Release Locking Ring (SH438)
9. Lever (SH881/G4)	23. Screw x 5 Per (SP989)
10. Slide Linkage Pin (SN220)	24. Left Handle (SH880L/G4)
11. Repeater Pawl (SH883)	25. 50ml Piston Seal Ring (SP346)
12. Pawl Return Spring (RP475)	26. 25ml Piston Seal Ring (SP887)
13. Lever Arm (SH882)	27. 50ml Push Rod Assembly (SXN229)
14. Return Spring (SP1075)	28. 25ml Push Rod Assembly (SXN385)
	29. Push Rod End Knob (SH1175)

PLEASE ORDER BY KIT AND PART NAME

SPARE PARTS AVAILABLE:

KIT No.	PART NAME	PART No.
WX1111	50ml Cylinder Assembly	1,2,3,4,7,25.
*WX1719	50ml Cylinder Assembly Prop	Not illustrated.
WX1433	50ml O-rings	25 (pack of five).
WX999	25ml Cylinder Assembly	1,2,3,4,7,25.
*WX1718	25ml Cylinder Assembly Prop	Not illustrated.
WX1434	25ml O-rings	26 (pack of five).
WX1421	4-Pack Coloured Knobs	Picture above.



njphillips.com

Our expertise is in your hands.

NJ PHILLIPS PTY LIMITED ABN 36 000 082 002
44 Gindurra Road, Somersby NSW, Australia 2250
Postal Address: Locked Bag 8, Gosford NSW, Australia 2250
Phone: +61 2 4340 2044 **Fax:** +61 2 4340 1991
Email: njp1@njphillips.com.au



The information contained in this document is legally privileged and confidential intended only for the use of the individual or entity. If the receiver of this document is not the intended recipient, the receiver is hereby notified that any dissemination, distribution or copy of this information is strictly prohibited. If this information is received in error, please notify us immediately. As our policy is one of continuous improvement, the manufacturer reserves the right to alter these specifications at any time. All products produced by NJ Phillips Pty Limited are identified by a unique batch number which is allocated to the product to allow traceability by the manufacturer and must not be removed if product integrity is to be maintained. Not all part numbers may be available.